

SUPPLEMENTARY DATA

Supplementary Table 1. The primers used in real-time polymerase chain reactions.

Col I α 1	sense	GTA TCA CCA AAC TCA GAA G
	antisense	ATC ATA GCC ATA GGA CAT C
Col I α 2	sense	CCA GGA AGT CCA GGA GTC
	antisense	TGA CAA AGG AGT TTC ATC TG
Col III α 1	sense	GGA GGG AGA AGT GGA GAG GAA AGG
	antisense	ATG GCT TCA AGG TGG CGT TAT CAG
Col IV α 1	sense	AAA GGG TCG GTG GGA GAG
	antisense	CGG GAA AGC CAG GGT TAC
GAPDH	sense	CGT GCC GCC TGG AGA AAC
	antisense	AGT GGG AGT TGC TGT TGA AGT

The primers for F4/80 and TNF α were purchased from Applied Biosystems.