

SUPPLEMENTARY DATA

Supplementary Table 1. Sequence of primers used for real-time PCR

Gene	Sequence of PCR primers
FoxP3	Sense: 5'-CACCTGGCTGGGAAAATGG-3'
	Antisense: 5'-GGAGCCCTTGTCGGATGAT-3'
CTLA-4	Sense: 5'-CCGAACTAACTGCTGCAAGGA-3'
	Antisense: 5'-CCCAGATTTATGTAATTGATCCAGAA-3'
GITR	Sense: 5'-GTTTTGGCTTCCAGTGTATCGA-3'
	Antisense: 5'-AACACAGTGAGAAACCCGAACT-3'
CD62L	Sense: 5'-CTCCTTGCCAGCCAAATGATAA-3'
	Antisense: 5'-CCTCTTCATTCCAGTGGCAGTC-3'
TGF- β	Sense: 5'-TGGAACCCACAACGAAATC-3'
	Antisense: 5'-GGGTTCAAGTACCGCTTCTC-3'
IL10	Sense: 5'-GCC TAA CAT GCT TCG AGA TC-3'
	Antisense: 5'-TGA TGT CTG GGT CTT GGT TC-3'
IL2	Sense: 5'-GAACTAAAGGGATCTGAAACAACATTC-3'
	Antisense: 5'-TGTTGAGATGATGCTTTGACAAAA-3'
IL4	Sense: 5'-AGAAGACTCTGTGCACCGAGTTGA-3'
	Antisense: 5'-CTCTCATGATCGTCTTTAGCCTTT-3'
IL17	Sense: 5'-ACTGCTACTGCTGCTGAGCCT-3'
	Antisense: 5'-TGGGCATCCTGGATTTCG-3'
IFN- γ	Sense: 5'-GGCATTTTGAAGAATTGGAAAG-3'
	Antisense: 5'-TTTGATGCTCTGGTCATCTT-3'
Granzyme A	Sense: 5'-GAGACTCGTGCAATGGAGATTCT-3'
	Antisense: 5'-AAGTGACCCCTCGGAAAACA-3'
GAPDH	Sense: 5'-TGCACCACCAACTGCTTAGC-3'
	Antisense: 5'-GGCATGGACTGTGGTCATGAG-3'